

HSC Safety Tips No. (20)

- **1. ALWAYS AVOID** potential burns or cuts while handling or sorting hot sterilized items or sharp instruments, and use caution when removing them from autoclaves/sterilizers or from steam lines that service the autoclaves.
- **2. ALWAYS** wear appropriate PPE when working with autoclaves/sterilizers. Using oven mitts for handling hot items, and steel-mesh gloves for handling or sorting sharp instruments are examples of appropriate PPE.
- 3. The autoclave/sterilizer door MUST ALWAYS be closed and locked before beginning the cycle.
- **4. DO NOT** remove items from an autoclave/sterilizer until they have cooled.
- **5. ALWAYS** use forceps or other tools to remove sharp instruments from baskets and autoclaves.

